**PRINT SPECIFICATION PAGE ONLY.
DO NOT PRINT THIS PAGE IN FINAL DFU.**

Pentron Clinical Marketing

Account: Pentron
Account Name: Pentron
Product Name: Simile

Part Number: DMC1657_Rev0

Description: Simile DFU

Physical specifications:

- a. Binding: Fold, stitch
- b. Corner Radius: N/A
- c. Colors: Black
- d. Size: Flat: 4.875" x 5.75"
Folded: 2.4375" x 5.75"
- e. Material: 20# bond white
- f. Adhesive material: N/A
- g. Coating: N/A
- h. Unwind direction: N/A

Revision Number: 0

Revision Date: 09/2014

Revision Summary: EC Rep Project

Languages

English	French	German
Italian	Spanish	Czech

Simile™

Universal Composite

English

INSTRUCTIONS FOR USE:

English

INSTRUCTIONS FOR USE

Français

French

MODE D'EMPLOI

Deutsch

German

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Italiano

Italian

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Español

Spanish

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čeština

Czech

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Indications:

- Simile Nano Hybrid Composite is a dental restorative material intended for use to restore carious lesions or structural defects in teeth in combination of conditioners such as bonding, luting, etching agents commonly used in tooth restoration.

Recommended for Class I, II, III, IV, V and VI Restorations

General Precautions

Avoid using the following:

- Hemodent with ferric sulfate or aluminum chloride.
- Retraction cord impregnated with the same.
- Pumice which contains oil.
- Compressed air which contains moisture or oil.

Preparation*

- Isolate with rubber dam or other suitable method.
 - Pumice tooth as needed.
 - Prepare the restorative site.
 - Place mylar strip(s) and wedge(s) where appropriate
- *In Class II Preparations, place a 0.5 -1mm enamel bevel along the cavosurface margins of the facial and lingual vertical walls of the preparation to improve enamel adhesion.*

Applying Bonding Agent

Please refer to the instructions of the specific bonding agent you are using.

Applying Simile nano-hybrid composite

- **For lighter shades** (A1, A2, B1, B2 and Incisal), place and cure in incremental layers from 0.5mm up to 3mm.
- **For darker shades** (A3, A3.5, B3, C2 and Universal Opaque), place and cure in incremental layers from 0.5mm up to 2mm.
- Pat and condense each layer with a smooth instrument using light pressure. For standard curing lights, cure each layer for 10 seconds. Cure final restoration an additional 20 seconds per surface. For high output halogen lights, cure each layer for 3 seconds and cure final restoration an additional 6 seconds per surface.
- Remove mylar strip, matrix band or crown form and check the restoration.

Trimming, Finishing, & Polishing

- Gross contouring may be accomplished using high speed diamond and 12 or 16 fluted carbide burs.
- Finer contouring/surface finishing may be accomplished with coarse Fini™ polishing disks, rubber points/ cups and finishing strips to create a smoother surface. To achieve a final polish, use the Fini polishing disks, starting with the coarse grits and moving down to the super fine grit. Other commercially available polishing systems can also be used.

Simile Nano-Hybrid Composite is offered in syringe and single dose delivery. Single dose delivery is intended for single patient (single use) only to prevent cross-contamination between patients.

NOTE: Refer to outer package for expiration date. Do not use after expiration date.

CAUTION: Uncured methacrylate resin may cause contact dermatitis and damage the pulp. Avoid contact with skin, eyes and soft tissue. Wash thoroughly with water after contact.

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Store at ambient temperature.

LIMITED WARRANTY —

Pentron's technical advice, whether verbal or in writing, is designed to assist dentists in using Pentron's product. The dentist assumes all risk and liability for damages arising out of the improper use of Pentron's product.

In the event of a defect in material or workmanship, Pentron will replace the defective product or part thereof, or reimburse the actual cost of the defective product. In order to take advantage of this limited warranty, the defective product must be returned to Pentron.

PENTRON

Manufactured by Pentron Clinical

1717 West Collins Ave. Orange, CA 92867 USA

800-551-0283 / 714-516-7557 / www.pentron.com

EC REP European Union Representative

SpofaDental a.s.

Markova 238

CZ-506 01 Jičín

Phone: +420 493 583 111